



Product Datasheet

Anti-TPM4 Antibody (orb1880615)


Description	Rabbit polyclonal antibody to TPM4
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-terminal region of human TPM4. The exact sequence is proprietary.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P67936
Entrez	7171
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/100), FC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of TPM4 expression...



Immunohistochemical analysis of TPM4 staining in tissue sections.



Flow cytometric analysis of U251 cells stained with Anti-TPM4 antibody.