


Product Datasheet

Lime1 antibody (orb18807)

Description	Goat polyclonal to Lime1
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Lime1 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 32kDa band observed in Mouse Heart lysates (predicted MW of 31kDa according to NP_076173.1). Recommended for use at 0.5-2 µg/ml.
Protein Sequence	ENVYESIKEMGL
Clonality	Polyclonal
MW	29.6
Dilution Range	ELISA: 1:128000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



orb18807
staining
(1ug/ml) of
Mouse Hear...