



Product Datasheet

Anti-TUBB2B Antibody (orb1880728)

Description	Rabbit polyclonal antibody to TUBB2B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-terminal region of human TUBB2B. The exact sequence is proprietary.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9CWF2 , Q3KRE8 , Q9BVA1
Entrez	347733 , 73710 , 291081
Dilution Range	WB (1/1000 - 1/2000), IH (1/50 - 1/200)
Expiration Date	12 months from date of receipt.



Western blot analysis of TUBB2B expressi...



Immunohistochemical analysis of TUBB2B
S...