



## Product Datasheet

### Anti-DCAF8L2 Antibody (orb1880739)

|                            |                                                                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to DCAF8L2                                                                                                      |
| <b>Reactivity</b>          | Human                                                                                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                                                                                               |
| <b>Tested Applications</b> | IF, IH, WB                                                                                                                                 |
| <b>Immunogen</b>           | KLH-conjugated synthetic peptide encompassing a sequence within the N-terminal region of human DCAF8L2. The exact sequence is proprietary. |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                  |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                           |
| <b>Note</b>                | For research use only                                                                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                                                                                 |
| <b>Source</b>              | Rabbit                                                                                                                                     |
| <b>Uniprot ID</b>          | <b>P0C7V8</b>                                                                                                                              |
| <b>Entrez</b>              | <b>347442</b>                                                                                                                              |
| <b>Dilution Range</b>      | WB (1/500 - 1/1000), IH (1/50 - 1/200), IF/IC (1/10 - 1/50)                                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                            |




Western blot analysis of DCAF8L2 express...



Immunohistochemical analysis of DCAF8L2

...



Immunofluorescent analysis of Anti-DCAF8...