



Product Datasheet

PDLIM1 antibody (orb18809)

Description	Goat polyclonal antibody to PDLIM1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CLIM1 / PDLIM1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 37kDa band observed in human lung and heart lysates (calculated MW of 36.1kDa according to NP_066272.1). Recommended for use at 0.1-0.3 µg/ml.
Protein Sequence	TPPEGYEVVTVFPK
Clonality	Polyclonal
MW	36.1
Entrez	9124
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb18809
staining
(0.1µg/ml) of
human l...



orb18809
(5µg/ml)
staining of
paraffin ...