


Product Datasheet

Human FGFR2 alpha ELISA Kit (orb1881034)

Description	Sandwich ELISA for quantitative detection of Human FGFR2 alpha concentrations.
Reactivity	Human
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P21802
Sensitivity	80 pg/ml
Entrez	2263
Expiration Date	Please enquire.



Representative
Standard
Curve