


Product Datasheet

Phosphatidylcholine Microplate Assay Kit (orb1881038)

Description	Detection and Quantification of Phosphatidylcholine Content.
Reactivity	All
Range	0.2 mmol/L - 20 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Sample Types	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	8002-43-5
Expiration Date	Please enquire.



Representative
Standard
Curve