


Product Datasheet

D-Amino Acid Microplate Assay Kit (orb1881039)

Description	Detection and Quantification of D-Amino Acid Content.
Reactivity	All
Range	0.5 mmol/L - 50 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Expiration Date	Please enquire.



Representative
Standard
Curve