



Product Datasheet

Butyryl Cholinesterase Microplate Assay Kit (orb1881040)

Description	Detection and Quantification of Butyryl Cholinesterase Activity.
Reactivity	All
Range	0.1 mmol/L - 5 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Sample Types	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9001-08-5
Expiration Date	Please enquire.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...