



Product Datasheet

RANBP1 antibody (orb18812)

Description	Goat polyclonal to RANBP1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	RANBP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 28-30kDa band observed in Hela and HepG2 lysates (predicted MW of 23kDa according to NP_002873). Recommended for use at 0.1-1 µg/ml.
Protein Sequence	LSVKEETKEDAEKQ
Clonality	Polyclonal
MW	23.3; 31.9; 23.2; 17.5
Entrez	5902
Dilution Range	ELISA: 1:16000, WB: 0.05-0.1 µg/ml, IHC-P: 3-5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HELA lysate using RANBP1 antibody.



Immunohistochemical staining of Human Kidney tissue using RANBP1 antibody.