



Product Datasheet

CD303 Antibody (PE) (orb1881234)


Description	CD303 is an approximately 38 kDa type II transmembrane glycoprotein with an extracellular C-type lec...
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	CD303 ectodomain fused with human IgG Fc domain
Target	CD303
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Uniprot ID	Q8WTT0
Entrez	170482
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD123
positive
CD303...



Flow
cytometry
surface
staining
pattern ...