


Product Datasheet

CD49c Antibody (APC) (orb1881239)

Description	CD49c / Integrin alpha 3 is a type I transmembrane glycoprotein proteolytically cleaved into two dis...
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	Human SSC-9 cell line (squamous cell carcinoma)
Target	CD49c
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Uniprot ID	P26006
Entrez	3675
Expiration Date	12 months from date of receipt.



Separation of HeLa cells stained using a...