


Product Datasheet

MAML1 antibody (orb18813)

Description	Goat polyclonal antibody to MAML1
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	MAML1 / Mastermind
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 110kDa band observed in Human Skeletal Muscle lysates (calculated MW of 108kDa according to NP_055572.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	VLPTCPMAEFALPR
Clonality	Polyclonal
MW	108
Entrez	9794
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 2.5-3.8µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Skeletal ...



Immunohistochemical staining of Human Pl...