

Product Datasheet

FOXP4 antibody (orb18814)

Description	Goat polyclonal antibody to FOXP4
Species/Host	Goat
Reactivity	Bovine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FOXP4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Customer's experiments gave a band at approx 75kDa in lysates of FOXP4-transfected HEK293 (calculated MW of 73.5kDa according to NP_001012426.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	VESASETIRSAPSG
Clonality	Polyclonal
MW	73.5; 72.2; 73.2
Entrez	116113
Dilution Range	ELISA: 1:4000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.