

Product Datasheet

Mouse anti-SARS-CoV2 (S-RBD) (XBB.1) antibody ELISA Kit (orb1881547)

Description	Mouse anti-SARS-CoV2 (S-RBD) (XBB.1) antibody ELISA Kit
Reactivity	Mouse
Range	3.125-200ng/ml
Target	anti-SARS-CoV2(S- RBD)(XBB.1)antibody
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	1.875ng/ml
Expiration Date	Please enquire.