


Product Datasheet

BRF2 antibody (orb18816)

Description	Goat polyclonal antibody to BRF2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	BRF2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 50kDa band observed in Hela and Human Placenta lysates (predicted MW of 48kDa according to NP_060780). Recommended for use at 1-3 µg/ml.
Protein Sequence	QAARQAATSVNPP
Clonality	Polyclonal
MW	46.5
Entrez	55290
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela lysate using...