

Product Datasheet

Recombinant Human Receptor-type tyrosine-protein phosphatase delta (PTPRD), partial (orb1881824)

Description

Recombinant Human Receptor-type tyrosine-protein phosphatase delta

Preservatives

If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Form/Appearance

Liquid or Lyophilized powder

Storage

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Tag

C-terminal 10xHis-tagged

Note

For research use only

Application notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Protein Sequence

ETPPRFTRTPVDQTVSGGVASFICQATGDPRPKIVWNKKGKKVSNQRFEVIEFD
DGSGSVLRIQPLRTPRDEAIYECVASNNVGEISVSTRLTVLREDQIPRGFPTIDMG
PQLKVVTRTATMLCAASGNPDPEITWFKDFLPVDTSNNGRIKQLRSESIGGT
PIRGALQIEQSEESDQGYECVATNSAGTRYAPANLYVRELREVRVPPRFSIPP
TNHEIMPGGSVNITCVAVGSPMPYVKWMLGAEDLTPEDDMPIGRNVLELNDVR
QSANYTCVAMSTLGVIEAIAQITVKALPKPPGTPVVTTESTATSITLTWDSGNPEPV
SYIIQHKKPKNSEELYKEIDGVATTRYVAGLSPYSDYEFVAVVNNIGRGPPSEP
VLTQTSEQAPSSAPRDVQARMLSTTILVQWKEPEEPNGQIQGYRVYTTMDPTQ
HVNWMKHNVADSQITTIGNLVPQKTSVKVLAFTSIGDGPLSSDIQVITQTGVP
GQPLNFKAEPESETSILLSWTPPRSDTIANVELVYKDGEHGEEQRITIEPGTSYRLQ
GLKPNLSLYFRLAARSPQGLGASTAEISARTMQSKPSAPPQDISCTSPSSSILVS
WQPPPVEKQNGIITEYSIKYTAVDGEDDKPHEILGIPSDTTKYLLEQLEKWTEYRIT
VTAHTDVGPPELSVLIRTNEDVPSGPPRKEVEAVNSTSVKVSWSRSPVPNKQ
HGQIRGYQVHYVRMENGEPKGQPMKDVMLADAQWFEFDDTTEHDMIISGLQP
ETSYSLTVTAYTTKGDGARSKPKLVSTTGAVPGKPRLVINHTQMNTALIQWHPPV
DTFGPLQGYRLKFKRKDMEPLTTLEFSEKEDHFTATDIHKGASYVFRLSARNKVG
FGEEMVKEISIPPEVPTGFPQNLHSEGTSTSVQLSWQPPVLAERNGIITKYTLTYR
DINIPLLPMEQLIVPADTMTLTGLKPDTTYDVKVRAHTSKGPGPYSPSVQFRTLQ
VDQVFAKNFHVAVMKTSVLLSWEIPENYNSAMPFKILYDDGKMVEEVDGRATQ
KLIVNLKPEKSYFVLNTRGNSAGGLQHRVTAKTAPDVLRTKPAFIGKTNLGDMIT
VOLPEVPANENIKGYIIIVPLKKSARGKFIKPWESPDEMEDELLEISKRKRSIRYG
REVELKPYIAAHFDVLPTEFTLGDDKHGGFTNKQLQSGQEYVFFVLAVMEHAES
KMYATSPYSDPVVMDLDPQIPITDEEE

Purity

Greater than 90% as determined by SDS-PAGE.

Source

Mammalian cell

MW

139.5 kDa

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558