

## Product Datasheet

# Recombinant *Saccharomyces cerevisiae* Actin-related protein 5 (ARP5), partial (orb1881937)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Recombinant <i>Saccharomyces cerevisiae</i> Actin-related protein 5 (ARP5), partial                                                                                                                                                                                                                                                                                                                         |
| <b>Preservatives</b>     | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.                                                                                                                                                                          |
| <b>Form/Appearance</b>   | Liquid or Lyophilized powder                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                             |
| <b>Tag</b>               | N-terminal 6xHis-tagged                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Application notes</b> | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Protein Sequence</b>  | GAKRNRHQATIDNDPNDTFGANDEDWLIYTDITQNPEAFEEALEYKYKDIVERLLLEHDPNFTEEDTLEAQYDWRNSILHLFLRGPRPHDSENIHEQHQMHLNVERIRVPEVIFQPTMGGQDQAGICELSETILLKKFGSQPGKLSQTSIDMVNNVLITGGNAKVPGLKERIVKEFTGFLPTGTNITVMSSDPSLDAWKGMAALARNEEQYRKTVISKKEYEEYGPEYIKEHKLGNTRYFED                                                                                                                                                          |
| <b>Purity</b>            | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Source</b>            | <i>E. coli</i>                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>MW</b>                | 32.2                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Uniprot ID</b>        | <b>P53946</b>                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                                                                                                                                                                                                                                                                                                                              |