


Product Datasheet

KLF3 antibody (orb18820)


Description	Goat polyclonal antibody to KLF3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	KLF3 / BKLF
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 40kDa band observed in Human Heart lysates (predicted MW of 40kDa according to NP_057615). Recommended for use at 1-3 µg/ml. Please note that the 40kDa Band is sometimes accompanied by a 80-90kDa band of unknown identity.
Protein Sequence	LMFDPVPVKQEAMD
Clonality	Polyclonal
MW	38.8
Entrez	51274
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lys...



Immunofluorescence analysis of A431 cell...



Immunofluorescence analysis of A549 cell...