

Product Datasheet

Anti-Calnexin/CANX Antibody Picoband (monoclonal, 8D10B3) (orb1882106)

Description	Anti-Calnexin/CANX Antibody Picoband (monoclonal, 8D10B3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human Calnexin/CANX recombinant protein (Position: E68-R582).
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	,,,,
Expiration Date	12 months from date of receipt.