

Product Datasheet

Anti-SLC4A1/CD233/Band 3 Antibody Picoband (monoclonal, 5G2G7) (orb1882129)

Description	Anti-SLC4A1/CD233/Band 3 Antibody Picoband (monoclonal, 5G2G7)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	E.coli-derived human SLC4A1 (Position: E28-N365). Human SLC4A1 shares 75.7% and 74.5% amino acid (aa) sequence identity with and rat SLC4A1, respectively.
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.