

## Product Datasheet

### Anti-ALPP Antibody Picoband (monoclonal, 6C7G3) (orb1882164)

|                            |                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Anti-ALPP Antibody Picoband<br>(monoclonal, 6C7G3)                                                               |
| <b>Species/Host</b>        | Mouse                                                                                                            |
| <b>Reactivity</b>          | Human                                                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                                     |
| <b>Tested Applications</b> | FC, ICC, IF, IHC, WB                                                                                             |
| <b>Immunogen</b>           | A synthetic peptide corresponding to a sequence at the C-terminus of human ALPP.                                 |
| <b>Preservatives</b>       | Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na <sub>2</sub> HPO <sub>4</sub> , 0.05mg NaN <sub>3</sub> . |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                      |
| <b>Concentration</b>       | Adding 0.2 ml of distilled water =f 500 µg/ml.                                                                   |
| <b>Storage</b>             | -20°C                                                                                                            |
| <b>Note</b>                | For research use only                                                                                            |
| <b>Isotype</b>             | IgG2b                                                                                                            |
| <b>Clonality</b>           | Monoclonal                                                                                                       |
| <b>Dilution Range</b>      | ,,,,                                                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                  |