

Product Datasheet

Anti-GAD65/GAD2 Antibody Picoband (monoclonal, 7G2) (orb1882274)

Description	Anti-GAD65/GAD2 Antibody Picoband (monoclonal, 7G2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human GAD65, different from the related mouse and rat sequences by one amino acid.
Preservatives	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg Na ₃ .
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water =f 500 µg/ml.
Storage	-20°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	,,,,
Expiration Date	12 months from date of receipt.