

Product Datasheet

Mouse anti-Hu IgE Purified, clone BE5 (Monoclonal) (orb1882385)

Description	Mouse anti-Hu IgE Purified, clone BE5 (Monoclonal)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified human IgE.
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Form/Appearance	Purified by protein-A affinity chromatography.
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,
Expiration Date	12 months from date of receipt.