


Product Datasheet

Pleckstrin antibody (orb18824)

Description	Goat polyclonal to Pleckstrin
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	Pleckstrin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 45kDa band observed in Human Heart and Peripheral Blood Mononucleocyte lysates (calculated MW of 40.1kDa according to NP_002655.2). Recommended for use at 0.05-1 µg/ml.
Protein Sequence	EWIKAIQMASRTGK
Clonality	Polyclonal
MW	40.1
Entrez	5341
Dilution Range	ELISA: 1:16000, WB: 0.01-0.3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb18824
staining
(0.05ug/ml)
of Human P...