

Product Datasheet

Mouse anti-Hu CD79b Purified, clone CB3-1 (Monoclonal) (orb1882408)

Description	Mouse anti-Hu CD79b Purified, clone CB3-1 (Monoclonal)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	Fraction of Ig-associated molecules isolated from Ramos B cells
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Form/Appearance	Purified by protein-A affinity chromatography.
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1 kappa
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.