

Product Datasheet

Mouse anti-Hu CD66e Purified, clone CB30 (Monoclonal) (orb1882419)

Description	Mouse anti-Hu CD66e Purified, clone CB30 (Monoclonal)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IHC-P, IP
Immunogen	Human carcinoembryonic antigen (CEA; CEACAM5)
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Form/Appearance	Purified by protein-A affinity chromatography.
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.