

Product Datasheet

Mouse anti-Hu CD47 Purified, clone MEM-122 (Monoclonal) (orb1882437)

Description	Mouse anti-Hu CD47 Purified, clone MEM-122 (Monoclonal)
Species/Host	Mouse
Reactivity	Human, Porcine, Primate
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, WB
Immunogen	COS-7 (African green monkey) cells
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Form/Appearance	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.