

Product Datasheet

Mouse anti-Hu CD34 Purified, clone 4H11[APG] (Monoclonal) (orb1882454)

Description	Mouse anti-Hu CD34 Purified, clone 4H11[APG] (Monoclonal)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Permanent human cell line derived from peripheral leucocytes of a patient suffering from chronic myeloid leukaemia.
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Form/Appearance	Purified by protein-A affinity chromatography.
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.