

Product Datasheet

Mouse anti-Hu CD274 Purified, clone 29E.2A3 (Monoclonal) (orb1882474)

Description	Mouse anti-Hu CD274 Purified, clone 29E.2A3 (Monoclonal)
Species/Host	Mouse
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr
Immunogen	Full length human CD274
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Form/Appearance	Purified by protein-A affinity chromatography.
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG2b kappa
Clonality	Monoclonal
Dilution Range	,
Expiration Date	12 months from date of receipt.