

Product Datasheet

Mouse anti-Hu CD18 Purified, clone MEM-48 (Monoclonal) (orb1882511)

Description	Mouse anti-Hu CD18 Purified, clone MEM-48 (Monoclonal)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IHC-P, IP, WB
Immunogen	Leukocytes of a patient suffering from a LGL-type leukemia.
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Form/Appearance	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	
Expiration Date	12 months from date of receipt.