

Product Datasheet

Mouse anti-Blood Group Lewis a Purified, clone 7LE (Monoclonal) (orb1882608)

Description	Mouse anti-Blood Group Lewis a Purified, clone 7LE (Monoclonal)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Mucin isolated by size-exclusion chromatography from an ovarian cyst fluid (pure endocervical type according to the Fenoglio's classification).
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Form/Appearance	Purified by protein-A affinity chromatography.
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.