

Product Datasheet

Annexin A10 antibody (orb18838)

Description

Goat polyclonal to Annexin A10

Species/Host

Goat

Reactivity

Bovine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Target

Annexin A10 / Annexin 14

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.IHC: In paraffin embedded Human Pancreas shows staining of islet of Langerhans . Recommended concentration, 3-5ug/ml.WB: Preliminary experiments (at 1 µg/ml) gave an approx 45-48kDa band in Hela, A431 and Human Kidney lysates. Please note that currently we cannot find an explanation in the literature for the band we observe given the predicted size of approx. 37kDa according to NP_009124. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

FCGDYVQGTIFPAPN

Clonality

Polyclonal

MW

37.3

Entrez


11199

Dilution Range

ELISA: 1:8000, WB: 1ug/ml, IHC-P: 3-5µg/ml

Expiration Date

12 months from date of receipt.



orb18838
(3.8ug/ml)
staining of
paraffin...