


Product Datasheet

Pescadillo antibody (orb18841)

Description	Goat polyclonal to Pescadillo
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Pescadillo
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 70kDa band observed in HeLa lysate (calculated MW of 68kDa according to NP_055118). Recommended concentration: 1-3 µg/ml.
Protein Sequence	GGLEKKKYERGSA
Clonality	Polyclonal
MW	68
Entrez	23481
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 2.5-3.8µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using Pescadillo antibody.



Immunohistochemical staining of Human Testis using Pescadillo antibody.