



Product Datasheet

KLF1 antibody (orb18849)

Description	Goat polyclonal antibody to KLF1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	KLF1 / EKLF
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Protein Sequence	ATAETALPSISTLT
Clonality	Polyclonal
MW	38.2
Entrez	10661
Dilution Range	ELISA: 1:128000, IHC-P: 4-6 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Bo...



Immunohistochemical staining of Human Ce...