



Product Datasheet

MPP5 antibody (orb18856)


Description	Goat polyclonal to MPP5
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	PALS / MPP5
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Skeletal Muscle shows strong alternating staining of muscle fibres in longitudinal section. Recommended concentration, 5-10ug/ml.WB: Approx 75kDa band observed in Human Brain lysate (calculated MW of 77.3kDa according to NP_071919). Recommended concentration: 1-3 µg/ml.
Protein Sequence	NKASPPFPLISNAQ
Clonality	Polyclonal
MW	77.3; 73.4
Entrez	64398
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb18856 (1ug/ml)
staining of Human
Brai...



orb18856 (5ug/ml)
staining of paraffin
e...



Immunohistochemical
staining of Human
Sm...