

Product Datasheet

**HDAC4 + 5 + 9 antibody (FITC) [Out of stock]
(orb188594)**

Description

Rabbit polyclonal antibody against

Species/Host	Rabbit
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human HDAC4 + 5 + 9
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	119/122/111kDa
Uniprot ID	P56524, Q9UQL6, Q99N13, Q9UKV0, Q6NZM9, Q9Z2V6, Q99P99
Dilution Range	ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.