


Product Datasheet

SLC16A7 antibody (orb18866)

Description	Goat polyclonal antibody to SLC16A7
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	SLC16A7 / MCT2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 40-45kDa band observed in A549, lung epithelial (lung carcinoma) lysates (calculated MW of 52.2 kDa according to NP_004722). Recommended concentration: 1-3 µg/ml.
Protein Sequence	SNAQSVTSERETNI
Clonality	Polyclonal
MW	52.2
Entrez	9194
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 lysate using SLC16A7 antibody.



Immunohistochemical staining of Human Testis tissue using SLC16A7 antibody.