

Product Datasheet

Insulin Receptor (phospho-Tyr999) antibody (FITC) (orb188680)

Description

Insulin Receptor (phospho-Tyr999)

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Insulin Receptor around the phosphorylation site of Tyr999
Target	Phospho-Insulin Receptor (Tyr999)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	affinity purified by Protein A
MW	70kDa
CAS Number	3
Uniprot ID	P06213
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.