



## Product Datasheet

### CRKL antibody (orb18870)

|                            |                                                                                                                                                                                                                          |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Goat polyclonal antibody to CRKL                                                                                                                                                                                         |
| <b>Species/Host</b>        | Goat                                                                                                                                                                                                                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                             |
| <b>Tested Applications</b> | ELISA, FC, IF, IHC, WB                                                                                                                                                                                                   |
| <b>Target</b>              | CRKL                                                                                                                                                                                                                     |
| <b>Preservatives</b>       | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.                                                              |
| <b>Storage</b>             | Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                                                                                               |
| <b>Note</b>                | For research use only                                                                                                                                                                                                    |
| <b>Application notes</b>   | ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 37kDa band observed in K562 and Jurkat lysates (calculated MW of 33.8kDa according to NP_005198). Recommended concentration: 0.01-0.03 µg/ml. |
| <b>Protein Sequence</b>    | KIFDPQNPDENE                                                                                                                                                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                               |
| <b>MW</b>                  | 33.8                                                                                                                                                                                                                     |
| <b>Entrez</b>              | <a href="#">1399</a>                                                                                                                                                                                                     |
| <b>Dilution Range</b>      | ELISA: 1:16000, WB: 0.01-0.03 µg/ml, IHC-P: 2.5-3.8 µg/ml                                                                                                                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                          |



Immunohistochemical staining of Human Th...



Western blot analysis of K562 lysate usi...