


Product Datasheet

DKK2 antibody (orb18872)

Description	Goat polyclonal antibody to DKK2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	DKK2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 28kDa band observed in Human Heart lysate (calculated MW of 28.5kDa according to NP_055236). Recommended concentration: 1-3 µg/ml
Protein Sequence	ATYSSKARLHVCQKI
Clonality	Polyclonal
MW	28.4
Entrez	27123
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb18872
(2ug/ml)
staining of
Human
Hear...