



Product Datasheet

TAOK3 antibody (orb18876)

Description	Goat polyclonal antibody to TAOK3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	JIK / TAOK3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Prostate shows staining of cytoplasm of secretory cells in the glands. Recommended concentration, 2-3ug/ml.WB: Approx 100kDa band observed in human Liver lysate (calculated MW of 105kDa according to NP_057365). Recommended concentration: 0.3-1.0 µg/ml.
Protein Sequence	CHKKDHFVIRDEAGHGD
Clonality	Polyclonal
MW	105
Entrez	51347
Dilution Range	ELISA: 1:32000, WB: 0.3-1.0 µg/ml, IHC-P: 2-3µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human liver lys...



Immunohistochemical staining of Human Pr...