


Product Datasheet

SAMSN1 antibody (orb18878)

Description	Goat polyclonal antibody to SAMSN1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SAMSN1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 48kDa band observed in Jurkat Tcell lysate (calculated MW of 41.7kDa according to NP_071419.2). The observed size is consistent with observations with competing products. Recommended concentration: 0.5-1.5 µg/ml
Protein Sequence	IRLDTLSASLGRSC
Clonality	Polyclonal
MW	41.7
Entrez	64092
Dilution Range	ELISA: 1:64000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat lysate u...