


Product Datasheet

VAV3 antibody (orb18879)


Description	Goat polyclonal antibody to VAV3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	VAV3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 80kDa band observed in human bone marrow lysates (calculated MW of 97.8kDa according to NP_006104). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	C-SGEQGTLKLPEK
Clonality	Polyclonal
MW	97.8; 32.6
Entrez	10451
Dilution Range	ELISA: 1:32000, WB: 0.5-1.5 µg/ml, IHC-P: 7.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Bone Marrow lysates



Immunohistochemical staining of Human skin tissue



Immunohistochemical staining of Human spleen tissue