

Product Datasheet

PITPNM1 antibody (orb18880)

Description	Goat polyclonal antibody to PITPNM1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	PITPNM / PITPNM1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Customer's experiments gave a band at approx 160kDa in lysates of a Mouse Macrophage cell line after 2 µg/ml antibody staining (calculated MW of 135kDa according to NP_004901.1). In ICC the customer found staining of vesicular cytoplasm in the same cell line.
Protein Sequence	KARSISLKLDSEE
Clonality	Polyclonal
MW	135; 135
Entrez	9600
Dilution Range	ELISA: 1:8000, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.