


Product Datasheet

TIRAP antibody (orb18890)

Description	Goat polyclonal antibody to TIRAP
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TIRAP / Mal (Isoform b)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 30kDa band observed in Human Liver lysates (calculated MW of 25.3kDa according to NP_683708.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	EGEGERDSATVSDL
Clonality	Polyclonal
MW	25.3
Entrez	114609
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb18890
(1ug/ml)
staining of
Human Live...