



Product Datasheet

FTCD antibody (orb18891)

Description	Goat polyclonal antibody to FTCD
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	58KGolgi protein(Internal)/FTCD
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 60kDa band observed in Human Liver lysates (calculated MW of 58.9kDa according to NP_006648.1and NP_996848.1). Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	CLREQGRGKDQPGRL
Clonality	Polyclonal
MW	58.9; 58.9; 61.3
Entrez	10841
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...



Immunohistochemical staining of Human Li...