



Product Datasheet

FHL1 antibody (orb18892)

Description	Goat polyclonal antibody to FHL1
Species/Host	Goat
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide with sequence CNKRFVFHNEQVY, from the internal region of the protein sequence according to NP_001070830.1; NP_001274729.1.
Target	Fhl1/ Slim1 (mouse)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Testis shows fibrous staining in the cytoplasm of developing sperm cells.. Recommended concentration, 2-2ug/ml.WB: Approx 30kDa band observed in Human Muscle lysates (calculated MW of 31.9kDa according to NP_001440). Recommended concentration: 0.05-0.3 µg/ml.
Protein Sequence	CNKRFVFHNEQVY
Clonality	Polyclonal
MW	31.9; 33.6
Entrez	2273
Dilution Range	ELISA: 1:32000, WB: 0.05-0.3 µg/ml, IHC-P: 2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Muscle ly...



Immunohistochemical staining of Human Te...