



Product Datasheet

EHD2 antibody (orb18895)


Description	Goat polyclonal antibody to EHD2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	EHD2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 60kDa band observed in Human Heart lysates (calculated MW of 61.2kDa according to NP_055416.2). In transfected HEK293 transiently expressing EHD2 a band of approx. 62kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 0.003-0.01 µg/ml.
Protein Sequence	CRLVPPSKRRHKGSA
Clonality	Polyclonal
MW	61.2
Entrez	30846
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



orb18895 Flow
cytometric analysis
of par...



orb18895
Immunofluorescence
analysis of ...



orb18895
Immunofluorescence
analysis of ...