



Product Datasheet

FBNP1 antibody (orb18896)

Description	Goat polyclonal antibody to FBNP1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	FBP17 / FBNP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 70kDa band observed in human testis lysates (calculated MW of 71.3kDa according to NP_055848.1). Recommended concentration: 0.25-2 µg/ml.
Protein Sequence	KQLESSKRRFERDC
Clonality	Polyclonal
MW	71.3
Entrez	23048
Dilution Range	ELISA: 1:32000, WB: 0.25-2 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human To...



Immunohistochemical staining of Human Sm...